


Mass Spectrum



Sample Information

Month-Day Processed : 03/23/23
Time Processed : 8:50:22 PM
Injection Volume : 0.1
Sample Name : Exendin-4
Cat No. : RP10874
Lot No. : E6347523B

Theoretical MW: 4186.66
Observed MW : 4187.1

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : ZJ22010150
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH