


Sample Name :GRGDSPC
Cat No. :RP20283
Lot No. :E6350726E
Time Processed :10:16:58
Month-Day-Year Processed :05/10/2026
Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	0
5.00	Pumps	B.Conc	0
20.00	Pumps	B.Conc	30
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	0
35.00	Pumps	B.Conc	0
35.01	Controller	Stop	

<<Column Performance>>
<Detector A>
Column :Inertsil ODS-SP 4.6 x 250 mm
Equipment: ZJ20010140

<Chromatogram>



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	10.117	4211100	128301	95.005
2	10.908	123168	15410	2.779
3	11.667	17079	1243	0.385
4	11.833	4504	910	0.102
5	12.800	76653	13702	1.729
Total		4432504	159566	100.000