

Sample Name: Glucagon-Like Peptide (GLP) II, human
 Cat No. : RP10769
 Lot No. : E6407322I
 Time Processed : 15:29:28
 Month-Day-Year Processed : 10/03/2022

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow: 1 ml/min
 Wavelength: 220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	Pump B Conc.	5
25.00	Pumps	Pump B Conc.	65
25.01	Pumps	Pump B Conc.	95
27.00	Pumps	Pump B Conc.	95
27.01	Pumps	Pump B Conc.	5
35.00	Pumps	Pump B Conc.	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column : Inertsil ODS-SP 4.6 x 250 mm
 Equipment: GK12010012

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	16.313	8525	865	0.105
2	17.117	7908	880	0.098
3	17.517	83365	1948	1.032
4	18.267	20989	1323	0.260
5	18.496	12246	1211	0.152
6	18.811	25657	1374	0.317
7	19.161	96962	7320	1.200
8	19.827	7763587	707944	96.063
9	20.383	25691	2264	0.318
10	20.608	17874	1406	0.221
11	22.922	18990	1140	0.235
Total		8081792	727676	100.000