


Mass Spectrum



Sample Information

|                       |                                       |                 |                 |
|-----------------------|---------------------------------------|-----------------|-----------------|
| Month-Day Processed : | 10/03/22                              | Equipment       | :ZJ22010150     |
| Time Processed :      | 3:24:38 PM                            | Interface Bias  | : +4.5 kV       |
| Injection Volume :    | 0.1                                   | Drying Gas Flow | :5 L/min        |
| Sample Name :         | Glucagon-Like Peptide (GLP) II, human | T.Flow          | :0.2 ml/min     |
| Cat No. :             | RP10769                               | B.conc          | :50%H2O/50%MeOH |
| Lot No. :             | E6407322I                             |                 |                 |
| Theoretical MW:       | 3766.20                               |                 |                 |
| Observed MW :         | 3766.8                                |                 |                 |