

Sample Name : Hepatitis B Virus Core (128-140)
Cat No. : RP17313
Lot No. : E6421522I
Time Processed : 11:04:15
Month-Day-Year Processed : 09/08/2022

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow: 1 ml/min
Wavelength: 220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm
Equipment: ZJ19010141

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	12.534	13351	1655	0.089
2	13.508	39691	4669	0.265
3	13.742	105556	11538	0.704
4	14.096	14702189	1177769	98.125
5	14.684	93700	7249	0.625
6	24.003	17403	874	0.116
7	24.513	11305	1001	0.075
Total		14983196	1204755	100.000