

Sample Name :Hepatitis B Virus Core (128-140)
 Cat No. :RP17313
 Lot No. :E6627322A
 Time Processed :1:29:28
 Month-Day-Year Processed :01/20/2022

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm
 Equipment: ZJ19010141

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	14.502	21041	1829	0.396
2	14.813	114844	16230	2.163
3	15.022	5091937	416766	95.892
4	15.868	4689	717	0.088
5	16.072	2795	426	0.053
6	16.283	22751	3785	0.428
7	16.404	35132	4512	0.662
8	16.717	16868	1909	0.318
Total		5310057	446173	100.000