


Mass Spectrum



Sample Information

Month-Day Processed :	01/19/22	Equipment	: GK11010007
Time Processed :	19:38:09	Interface Bias	: +4.5 kV
Injection Volume :	0.4	Drying Gas Flow	: 5 L/min
Sample Name :	Hepatitis B Virus Core (128-140)	T.Flow	: 0.2 ml/min
Cat No. :	RP17313	B.conc	: 50% H ₂ O / 50% MeOH
Lot No. :	E6627322A		
Theoretical MW :	1406.63		
Observed MW :	1406.4		