

Sample Name : Glucagon-Like Peptide (GLP) II, human

Cat No. : RP10769

Lot No. : E6630223H

Time Processed :21:40:28

Month-Day-Year Processed :08/07/2023

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.00	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: GR11010440

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	12.573	12736	1709	0.086
2	12.888	80578	23329	0.543
3	13.004	14672424	2683208	98.838
4	13.281	58685	9389	0.395
5	13.585	20524	4025	0.138
Total		14844946	2721660	100.000