


Mass Spectrum



Sample Information

|                       |                              |                     |           |                 |                 |
|-----------------------|------------------------------|---------------------|-----------|-----------------|-----------------|
| Month-Day Processed : | 08/15/26                     | Interface           | :ESI      | Equipment       | : ZJ21010035    |
| Time Processed :      | 9:32:32                      | Nebulizing Gas Flow | :1.5L/min | Interface Bias  | : +4.5 kV       |
| Injection Volume :    | 0.3                          | CDL Temp            | :250      | Drying Gas Flow | :5 L/min        |
| Sample Name :         | KRAS WT Peptide (VVVGAGGVGK) | Block Temp          | :200      | T.Flow          | :0.2 ml/min     |
| Cat No. :             | RP30235                      |                     |           | B.conc          | :50%H2O/50%MeOH |
| Lot No. :             | E6654326H                    |                     |           |                 |                 |
| Theoretical MW :      | 842.00                       |                     |           |                 |                 |
| Observed MW :         | 842.2                        |                     |           |                 |                 |