

Sample Name : Influenza A NP(366-374) Strain A/PR/8/35  
Cat No. :RP20267  
Lot No. : E6793922D  
Time Processed: 4:17:23  
Year-Month-Day: 2022/4/14

Pump A : 0.065% trifluoroacetic in 100% water (v/v)  
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)  
Total Flow:1 ml/min

<<Detector A>>  
Wavelength Ch1 : 220 nm

<<LC Time Program>>


| Time  | Module     | Command         | Value | Comment |
|-------|------------|-----------------|-------|---------|
| 0.01  | Pumps      | Solvent B Conc. | 5     |         |
| 25.00 | Pumps      | Solvent B Conc. | 65    |         |
| 25.01 | Pumps      | Solvent B Conc. | 95    |         |
| 27.00 | Pumps      | Solvent B Conc. | 95    |         |
| 27.01 | Pumps      | Solvent B Conc. | 5     |         |
| 33.00 | Pumps      | Solvent B Conc. | 5     |         |
| 33.01 | Controller | Stop            |       |         |

<<Column Performance>>

<Detector A>  
Column : Inertsil ODS-SP 4.6 x 250 mm  
Equipment:ZJ19010015

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

| Peak# | Ret. Time | Area    | Height | Area%   |
|-------|-----------|---------|--------|---------|
| 1     | 5.184     | 33112   | 3351   | 0.988   |
| 2     | 5.702     | 26889   | 2499   | 0.802   |
| 3     | 6.567     | 11306   | 813    | 0.337   |
| 4     | 6.990     | 3202828 | 332006 | 95.531  |
| 5     | 7.778     | 5971    | 711    | 0.178   |
| 6     | 8.011     | 11609   | 1414   | 0.346   |
| 7     | 12.897    | 16414   | 1798   | 0.490   |
| 8     | 13.810    | 44542   | 3193   | 1.329   |
| Total |           | 3352670 | 345785 | 100.000 |