


Mass Spectrum



Sample Information		Equipment	: GK11010007
Month-Day Processed :	04/14/22	Interface Bias	: +4.5 kV
Time Processed :	8:45:06	Drying Gas Flow	: 5 L/min
Injection Volume :	0.2	T.Flow	: 0.2 ml/min
Sample Name :	Influenza A NP(366-374) Strain A/PR/8/35	B.conc	: 50%H2O/50%M
Cat No. :	RP20267		
Lot No. :	E6793922D		
Theoretical MW :	1026.12		
Observed MW :	1025.8		