


Mass Spectrum



Sample Information

Month-Day Processed :	04/30/22	Equipment	: GK11010007
Time Processed :	18:29:21	Interface Bias	: +4.5 kV
Injection Volume :	0.2	Drying Gas Flow	: 5 L/min
Sample Name :	β-Amyloid (25-35)	T.Flow	: 0.2 ml/min
Cat No. :	RP10008	B.conc	: 50% H2O/50% MeOH
Lot No. :	E6795522D		
Theoretical MW :	1060.27		
Observed MW :	1060.4		