


Mass Spectrum



Sample Information
Month-Day Processed : 01/26/22
Time Processed : 16:46:45
Injection Volume : 0.2
Sample Name : TAT peptide
Cat No. : RP20256
Lot No. : E6829322A
Theoretical MW : 1621.91
Observed MW : 1622.0

Interface : ESI
Nebulizing Gas Flow : 1.5 L/min
CDL Temp : 250
Block Temp : 200
Equipment : ZJ21010035
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH