


# Mass Spectrum



## Sample Information

Month-Day Processed : 04/16/26  
Time Processed : 8:06:54  
Injection Volume : 0.5  
Sample Name : RGD  
Cat No. : RP20297  
Lot No. : E6863426D

## Interface

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

## Equipment

Equipment : ZJ21010035  
Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50%H<sub>2</sub>O/50%MeOH

Theoretical MW : 1025.06  
Observed MW : 1024.8