

Sample Name :Influenza A NP(366-374) Strain A/PR/8/35

Cat No. : RP20267

Lot No. : E6887823F

Time Processed :4:23:08

Month-Day-Year Processed :06/18/2023

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	Pump A B.Conc	5
25.00	Pumps	Pump A B.Conc	65
25.01	Pumps	Pump A B.Conc	95
27.00	Pumps	Pump A B.Conc	95
27.01	Pumps	Pump A B.Conc	5
35.00	Pumps	Pump A B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: SS-CM-0286

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	8.500	7113981	562657	98.678
2	8.900	40201	3574	0.558
3	10.292	40383	4988	0.560
4	10.633	14706	2148	0.204
Total		7209271	573366	100.000