

Sample Name :KRAS G12V Peptide (VVVGAVGVGK)
Cat No. : RP30229
Lot No. : E6914724E

Time Processed :22:18:33
Month-Day-Year Processed :05/26/2024

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm
Equipment: ZJ19010141

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	10.214	21127	4165	1.059
2	10.479	1904380	335559	95.459
3	10.704	23707	3234	1.188
4	12.058	16165	1497	0.810
5	12.188	18097	1989	0.907
6	16.986	11505	1546	0.577
Total		1994979	347989	100.000