


Mass Spectrum



Sample Information
 Month-Day Processed : 07/22/26
 Time Processed : 13:56:35
 Injection Volume : 0.4
 Sample Name : Exendin-4
 Cat No. : RP10874
 Lot No. : E6942626G
 Theoretical MW : 4186.66
 Observed MW : 4186.4

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200
 Equipment : ZJ21010035
 Interface Bias : +4.5 kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2 ml/min
 B.conc : 50%H₂O/50%MeOH