

Sample Name: C-Type Natriuretic Peptide (CNP) (1-22), human
 Cat No. : RP11110
 Lot No. : E6994423G
 Time Processed : 12:15:27
 Month-Day-Year Processed : 07/18/2023

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow: 1 ml/min
 Wavelength: 220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column : Inertsil ODS-SP 4.6 x 250 mm

Equipment: GR11010440

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	14.597	25386	3373	0.193
2	14.850	93121	15340	0.709
3	15.025	116295	10380	0.885
4	15.395	12564416	1542130	95.655
5	15.563	153001	43243	1.165
6	17.252	9284	1670	0.071
7	17.818	173582	20512	1.322
Total		13135085	1636647	100.000