


Mass Spectrum



Sample Information		Interface		Equipment	
Month-Day Processed :	04/25/26	Nebulizing Gas Flow :	1.5L/min	Interface Bias	:ZJ22010150
Time Processed :	11:05:10 AM	CDL Temp :	250	Drying Gas Flow	: +4.5 kV
Injection Volume :	0.5	Block Temp :	200	T.Flow	:5 L/min
Sample Name :	Glucagon-Like Peptide (GLP) II, human			B.conc	:0.2 ml/min
Cat No. :	RP10769				:50%H2O/50%MeOH
Lot No. :	E7030726D				
Theoretical MW :	3766.20				
Observed MW :	3765.9				