

Sample Name:Glucagon-Like Peptide (GLP) II, human
Cat No. : RP10769
Lot No. : E7042723H
Time Processed :12:45:05
Month-Day-Year Processed :09/19/2023

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>


Time	Module	Command	Value
0.01	Pumps	Pump A B.Conc	5
25.00	Pumps	Pump A B.Conc	65
25.01	Pumps	Pump A B.Conc	95
27.00	Pumps	Pump A B.Conc	95
27.01	Pumps	Pump A B.Conc	5
35.00	Pumps	Pump A B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>

<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm
Equipment: GR11010440

<Chromatogram>



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	13.719	183654	26381	1.319
2	14.184	156915	22372	1.127
3	14.646	139113	14189	0.999
4	14.938	13409188	1426890	96.289
5	15.247	37073	6306	0.266
Total		13925943	1496137	100.000