


Mass Spectrum



Sample Information		Equipment	
Month-Day Processed :	09/19/23	Interface	: ZJ21010035
Time Processed :	10:34:54	Interface Bias	: +4.5 kV
Injection Volume :	0.3	Drying Gas Flow	: 5 L/min
Sample Name :	Glucagon-Like Peptide (GLP) II, human	T.Flow	: 0.2 ml/min
Cat No. :	RP10769	B. conc	: 50% H2O/50% MeOH
Lot No. :	E7042723H		
Theoretical MW:	3766.20		
Observed MW :	3765.9		