

Sample Name :C-Type Natriuretic Peptide (CNP) (1-22), human
Cat No. :RP11110
Lot No. : E7079122K
Time Processed :8:57:32 AM
Month-Day-Year Processed :11/16/2022

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm
Equipment: ZJ19010141

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	13.389	9734	1361	0.433
2	13.695	9518	848	0.423
3	13.961	2207339	290713	98.088
4	14.317	23786	1975	1.057
Total		2250377	294896	100.000