

Sample Name :Glucagon-Like Peptide (GLP) II, human

Cat No. :RP10769

Lot No. :E7146825G

Time Processed :13:34:03

Month-Day-Year Processed :08/04/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: SS-CM-0310

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	14.750	62945	9144	0.823
2	15.058	37936	4038	0.496
3	15.583	6415	1111	0.084
4	15.708	9414	808	0.123
5	16.233	6702	1941	0.088
6	16.458	7515808	1066072	98.218
7	16.650	12981	5169	0.170
Total		7652201	1088283	100.000