

Sample Name :Tau peptide (275-305) (Repeat2 domain)
Cat No. :RP30887
Lot No. :E7162926B
Time Processed :13:50:36
Month-Day-Year Processed :02/21/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm
Equipment: ZJ20010264

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	11.108	16679	3795	0.315
2	11.242	68372	13814	1.290
3	11.417	5131988	699336	96.803
4	11.792	68131	9358	1.285
5	22.267	16315	1635	0.308
Total		5301484	727938	100.000