

Sample Name:KRAS G12V Peptide (VVVGAVGVGK)

Cat No. : RP30229

Lot No. : E7183126A

Time Processed :0:39:02

Month-Day-Year Processed :01/16/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: ZJ24010222

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	8.075	72954	11523	1.360
2	8.750	5245281	518107	97.756
3	9.150	33377	4206	0.622
4	16.583	14065	1743	0.262
Total		5365677	535579	100.000