


Mass Spectrum



Sample Information

Month-Day Processed :	01/15/26	Equipment	: ZJ21010035
Time Processed :	19:40:50	Interface Bias	: +4.5 kV
Injection Volume :	0.3	Drying Gas Flow	: 5 L/min
Sample Name :	KRAS G12V Peptide (VVVGAVGVGK)	T.Flow	: 0.2 ml/min
Cat No. :	RP30229	B.conc	: 50%H2O/50%MeOH
Lot No. :	E7183126A		
Theoretical MW:	884.08		
Observed MW :	883.8		