


# Mass Spectrum



|                          |           |                     |                 |
|--------------------------|-----------|---------------------|-----------------|
| Sample Information       |           | Equipment           |                 |
| Month-Day Processed :    | 08/03/24  | Interface           | :ZJ22010150     |
| Time Processed :         | 15:57:31  | Nebulizing Gas Flow | : +4.5 kV       |
| Injection Volume :       | 0.2       | CDL Temp            | :5 L/min        |
| Sample Name :            | RGD       | Block Temp          | :0.2 ml/min     |
| Cat No. :                | RP20297   | Drying Gas Flow     | :50%H2O/50%MeOH |
| Lot No. :                | E7224224H | T.Flow              |                 |
| Theoretical MW : 1025.06 |           | B.conc              |                 |
| Observed MW : 1024.8     |           |                     |                 |