

Sample Name :C-Type Natriuretic Peptide (CNP) (1-22), human

Cat No. :RP11110

Lot No. :E7277326G

Time Processed :1:20:36

Month-Day-Year Processed :08/02/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	Solvent B Conc.	5
25.00	Pumps	Solvent B Conc.	65
25.01	Pumps	Solvent B Conc.	95
27.00	Pumps	Solvent B Conc.	95
27.01	Pumps	Solvent B Conc.	5
33.00	Pumps	Solvent B Conc.	5
33.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ19010015

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	12.050	19581	2883	0.539
2	12.183	7354	1233	0.202
3	12.300	7401	1176	0.204
4	12.558	3584630	478393	98.592
5	13.683	16856	1797	0.464
Total		3635821	485481	100.000