


Mass Spectrum



Sample Information

Month-Day Processed : 08/01/26
Time Processed : 19:44:40
Injection Volume : 0.3
Sample Name : C-Type Natriuretic Peptide (CNP) (1-22), human
Cat No. : RP11110
Lot No. : E7277326G
Theoretical MW : 2197.60
Observed MW : 2197.8

08/01/26
19:44:40
0.3
C-Type Natriuretic Peptide (CNP) (1-22), human

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : ZJ21010035
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H₂O/50%MeOH