

Sample Name :Influenza A NP(366-374) Strain A/PR/8/35
Cat No. : RP20267
Lot No. : E7331923K
Time Processed :1:07:19
Month-Day-Year Processed :11/17/2023

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ21010026

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	6.258	4089	3569	0.106
2	6.392	3825654	451737	99.168
3	6.633	6966	2616	0.181
4	7.542	21052	3160	0.546
Total		3857761	461081	100.000