

Sample Name:C-Type Natriuretic Peptide (CNP) (1-22), human
Cat No. : RP11110
Lot No. : E7437726F
Time Processed :2:31:36
Month-Day-Year Processed :06/23/2026


Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>
<Detector A>
Column :Inertsil ODS-SP 4.6 x 250 mm
Equipment: ZJ25010054

<Chromatogram>



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	12.875	15453	1953	0.447
2	13.050	10410	1440	0.301
3	13.417	3419325	413714	98.935
4	13.800	10955	1196	0.317
Total		3456143	418303	100.000