


Mass Spectrum



Sample Information	Interface	Equipment
Month-Day Processed : 06/22/26	Nebulizing Gas Flow : 1.5L/min	: ZJ22010150
Time Processed : 15:37:52	CDL Temp : 250	: +4.5 kV
Injection Volume : 0.2	Block Temp : 200	: 5 L/min
Sample Name : C-Type Natriuretic Peptide (CNP) (1-22), human		: 0.2 ml/min
		: 50% H2O/50% MeOH

Cat No. : RP11110
 Lot No. : E7437726F
 Theoretical MW: 2197.60
 Observed MW : 2197.2