

Sample Name :Hepatitis B Virus Core (128-140)
 Cat No. : RP17313
 Lot No. : E7478722G
 Time Processed :12:01:54 AM
 Month-Day-Year Processed :07/30/2022

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ21010026

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	13.018	8070	844	0.156
2	13.312	11219	1122	0.217
3	13.548	43025	4426	0.831
4	13.981	5094603	560280	98.386
5	24.735	21264	2299	0.411
Total		5178180	568971	100.000