

Sample Name :KRAS G12V Peptide (VVVGAVGVGK)

Cat No. : RP30229

Lot No. : E7495426C

Time Processed :9:59:46 AM

Month-Day-Year Processed :04/02/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ17010508

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	10.275	17186	2419	0.665
2	10.592	5812	1066	0.225
3	10.683	6414	1167	0.248
4	10.875	1796	404	0.069
5	11.033	2450	727	0.095
6	11.142	15963	1909	0.618
7	11.508	2456886	300704	95.057
8	11.883	25212	3446	0.975
9	12.708	52918	6483	2.047
Total		2584639	318325	100.000