

Sample Name:Melittin,honey bee
Cat No. : RP20415
Lot No. : E7516025F
Time Processed :9:16:09
Month-Day-Year Processed :06/25/2025


Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>
<Detector A>
Column :Inertsil ODS-SP 4.6 x 250 mm
Equipment: ZJ19010140

<Chromatogram>



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	19.642	21169	1509	0.354
2	19.900	57751	3314	0.965
3	20.408	5716129	491079	95.562
4	20.692	123621	13300	2.067
5	21.008	62903	4427	1.052
Total		5981574	513629	100.000