


Mass Spectrum



Sample Information

Month-Day Processed :	06/07/24	Equipment	: ZI21010035
Time Processed :	9:47:19	Interface Bias	: +4.5 kV
Injection Volume :	0.2	Drying Gas Flow	: 5 L/min
Sample Name :	β-Amyloid (25-35)	T.Flow	: 0.2 ml/min
Cat No. :	RP10008	B.conc	: 50%H2O/50%MeOH
Lot No. :	E7517824F		
Theoretical MW:	1060.27		
Observed MW :	1060.0		