

Sample Name:TAT peptide
 Cat No. : RP20256
 Lot No. : E7587422K
 Time Processed : 19:31:08
 Month-Day-Year Processed : 11/28/2022

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	Pump B Conc.	0
5.00	Pumps	Pump B Conc.	0
20.00	Pumps	Pump B Conc.	30
25.00	Pumps	Pump B Conc.	65
25.01	Pumps	Pump B Conc.	95
31.00	Pumps	Pump B Conc.	95
31.01	Pumps	Pump B Conc.	0
40.00	Pumps	Pump B Conc.	0
40.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment:GK12010012

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	10.921	176823	3759	2.081
2	11.574	8111658	854289	95.457
3	11.879	209204	9289	2.462
Total		8497685	867337	100.000