


Mass Spectrum



Sample Information		Interface		Equipment	
Month-Day Processed :	10/19/22	Nebulizing Gas Flow :1.5L/min		Interface Bias	
Time Processed :	20:42:45	CDL Temp :250		Drying Gas Flow	
Injection Volume :	0.1	Block Temp :200		T.Flow	
Sample Name :	DYKDDDDK Peptide			B.conc	
Cat No. :	RP10586				
Lot No. :	E7612322J				
Theoretical MW:1012.97					
Observed MW : 1012.8					