

Sample Name:Influenza A NP(366-374) Strain A/PR/8/35
Cat No. : RP20267
Lot No. : E7673326A
Time Processed :0:27:48
Month-Day-Year Processed :01/21/2026


Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>
<Detector A>
Column :Inertsil ODS-SP 4.6 x 250 mm
Equipment: ZJ20010140

<Chromatogram>



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	6.825	7828417	868568	96.454
2	7.100	52869	11424	0.651
3	7.583	20570	2997	0.253
4	7.725	9969	2610	0.123
5	7.842	86028	12062	1.060
6	15.400	118362	9965	1.458
Total		8116215	907626	100.000