

Sample Name: Tau peptide (275-305) (Repeat2 domain)
Cat No. : RP30887
Lot No. : E7807525H
Time Processed : 4:47:29
Month-Day-Year Processed : 08/03/2025


Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow: 1 ml/min
Wavelength: 220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>
<Detector A>
Column : Inertsil ODS-3 4.6 x 250 mm
Equipment: ZJ19010142

<Chromatogram>



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	7.567	26197	3187	0.178
2	9.275	59709	7979	0.406
3	9.458	258392	54031	1.759
4	9.708	14061035	1300024	95.696
5	10.092	44676	7054	0.304
6	10.200	35651	5653	0.243
7	10.600	75771	5232	0.516
8	11.342	131946	20359	0.898
Total		14693375	1403519	100.000