


Mass Spectrum



Sample Information

Month-Day Processed :	04/24/22	Equipment	: ZJ21010035
Time Processed :	7:36:51	Interface Bias	: +4.5 kV
Injection Volume :	0.2	Drying Gas Flow	: 5 L/min
Sample Name :	RGD	T.Flow	: 0.2 ml/min
Cat No. :	RP20297	B.conc	: 50% H_2O /50%MeOH
Lot No. :	E7879022D		
Theoretical MW :	1025.06		
Observed MW :	1025.0		