

Sample Name :C-Type Natriuretic Peptide (CNP) (1-22), human
 Cat No. : RP11110
 Lot No. : E7904425G
 Time Processed :4:17:12
 Month-Day-Year Processed :07/24/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ17010506

<Chromatogram>



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	11.300	14352	2533	0.378
2	11.775	3768414	473024	99.362
3	12.133	9840	1557	0.259
Total		3792607	477114	100.000