


Mass Spectrum



Sample Information

Month-Day Processed :	07/23/25	Interface	:ESI	Equipment	:ZJ22010150
Time Processed :	16:44:34	Nebulizing Gas Flow	:1.5L/min	Interface Bias	: +4.5 kV
Injection Volume :	0.1	CDL Temp	:250	Drying Gas Flow	:5 L/min
Sample Name :	C-Type Natriuretic Peptide (CNP) (1-22), human	Block Temp	:200	T.Flow	:0.2 ml/min
Cat No. : RP11110				B.conc	:50%H ₂ O/50%MeOH
Lot No. : E7904425G					
Theoretical MW: 2197.60					
Observed MW : 2197.5					