

Sample Name :Melittin,honey bee
Cat No. : RP20415
Lot No. : E7908325G
Time Processed :3:41:45
Month-Day-Year Processed :07/24/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ17010506

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	14.475	13278	1438	0.233
2	17.858	25643	2297	0.450
3	18.300	5469804	575781	96.017
4	18.650	126891	13043	2.227
5	18.842	61077	5269	1.072
Total		5696693	597828	100.000