

Sample Name:KRAS G12V Peptide (VVVGAVGVGK)
Cat No. : RP30229
Lot No. : E7930425G
Time Processed :21:51:39
Month-Day-Year Processed :07/25/2025


Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>
<Detector A>
Column :Inertsil ODS-SP 4.6 x 250 mm
Equipment: ZJ17010506

<Chromatogram>



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	6.817	37985	6077	1.726
2	7.042	9951	1287	0.452
3	7.317	6590	1527	0.299
4	7.492	2121759	253262	96.419
5	8.342	8688	1343	0.395
6	16.833	15580	1493	0.708
Total		2200553	264990	100.000