

Sample Name :Hepatitis B Virus Core (128-140)
 Cat No. : RP17313
 Lot No. : E7936325G
 Time Processed :1:49:49
 Month-Day-Year Processed :07/26/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: ZJ23010113

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	11.242	296726	41878	2.423
2	11.531	32497	2829	0.265
3	12.032	135789	11861	1.109
4	12.310	11759070	1540546	96.011
5	13.187	15487	1523	0.126
6	15.473	8093	927	0.066
Total		12247662	1599564	100.000