

Sample Name :Influenza A NP(366-374) Strain A/PR/8/35

Cat No. :RP20267

Lot No. : E7948725I

Time Processed :15:43:15

Month-Day-Year Processed :09/07/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

| Time  | Module     | Command | Value |
|-------|------------|---------|-------|
| 0.01  | Pumps      | B.Conc  | 5     |
| 25.00 | Pumps      | B.Conc  | 65    |
| 25.01 | Pumps      | B.Conc  | 95    |
| 27.00 | Pumps      | B.Conc  | 95    |
| 27.01 | Pumps      | B.Conc  | 5     |
| 35.00 | Pumps      | B.Conc  | 5     |
| 35.01 | Controller | Stop    |       |

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ21010376

### <Chromatogram>

mV



### <Peak Table>

Detector A Channel 1 220nm

| Peak# | Ret. Time | Area    | Height | Area%   |
|-------|-----------|---------|--------|---------|
| 1     | 5.992     | 75921   | 10049  | 2.269   |
| 2     | 6.492     | 3190775 | 421777 | 95.374  |
| 3     | 6.717     | 18576   | 3924   | 0.555   |
| 4     | 15.058    | 60266   | 6018   | 1.801   |
| Total |           | 3345538 | 441768 | 100.000 |