

Sample Name:KRAS G12V Peptide (VVVGAVGVGK)
Cat No. : RP30229
Lot No. : E7987724G
Time Processed :9:14:31
Month-Day-Year Processed :07/31/2024

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm
Equipment: ZJ23010112

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	10.000	6121	1725	0.360
2	10.258	1620442	309527	95.334
3	11.333	47299	6492	2.783
4	14.850	25896	3737	1.524
Total		1699758	321482	100.000