

Sample Name :Influenza A NP(366-374) Strain A/PR/8/35
 Cat No. :RP20267
 Lot No. :E7998522E
 Time Processed :3:02:10 PM
 Month-Day-Year Processed :05/27/2022

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	Pump A B.Conc	5
25.00	Pumps	Pump A B.Conc	65
25.01	Pumps	Pump A B.Conc	95
27.00	Pumps	Pump A B.Conc	95
27.01	Pumps	Pump A B.Conc	5
35.00	Pumps	Pump A B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm
 Equipment: SS-CM-0309

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	4.310	15061	2301	0.375
2	5.993	44040	6561	1.096
3	11.080	42619	6090	1.061
4	11.593	3850297	556885	95.837
5	11.779	65544	16849	1.631
Total		4017561	588687	100.000