


Mass Spectrum



Sample Information		Equipment : ZJ21010035	
Month-Day Processed :	04/11/23	Interface Bias :	+4.5 kV
Time Processed :	21:11:55	Drying Gas Flow :	5 L/min
Injection Volume :	0.2	T.Flow :	0.2 ml/min
Sample Name :	TAT peptide	B.conc :	50%H2O/50%MeOH
Cat No. :	RP20256		
Lot No. :	E8002323D		
Theoretical MW :	1621.91		
Observed MW :	1621.6		